## **Product data sheet**



MedKoo Cat#: 530478				
Name: Lavendustin B		он о		
CAS: 125697-91-8				
Chemical Formula: C <sub>21</sub> H <sub>19</sub> NO <sub>5</sub>		OH		
Exact Mass: 365.1263				
Molecular Weight: 365.385				
Product supplied as:	Powder			
Purity (by HPLC):	≥ 98%			
Shipping conditions	Ambient temperature			
Storage conditions:	Powder: -20°C 3 years; 4°C 2 years.	OH OH		
_	In solvent: -80°C 3 months; -20°C 2 weeks.			

#### 1. Product description:

Lavendustin B is a Tyrosine Kinase Inhibitor amd an inhibitor of HIV-1 integrase (IN) interaction with its cognate cellular cofactor, lens epithelium-derived growth factor (LEDGF/p75).

#### 2. CoA, QC data, SDS, and handling instruction

SDS and handling instruction, CoA with copies of QC data (NMR, HPLC and MS analytical spectra) can be downloaded from the product web page under "QC And Documents" section. Note: copies of analytical spectra may not be available if the product is being supplied by MedKoo partners. Whether the product was made by MedKoo or provided by its partners, the quality is 100% guaranteed.

3. Solubility data

Solvent	Max Conc. mg/mL	Max Conc. mM
DMF	10.0	27.37
DMSO	44.33	121.33
Ethanol	41.5	113.58

4. Stock solution preparation table:

ii brock bolation preparation table.					
Concentration / Solvent Volume / Mass	1 mg	5 mg	10 mg		
1 mM	2.74 mL	13.68 mL	27.37 mL		
5 mM	0.55 mL	2.74 mL	5.47 mL		
10 mM	0.27 mL	1.37 mL	2.74 mL		
50 mM	0.06 mL	0.27 mL	0.55 mL		

#### 5. Molarity Calculator, Reconstitution Calculator, Dilution Calculator

Please refer the product web page under section of "Calculator"

### 6. Recommended literature which reported protocols for in vitro and in vivo study

In vitro study

TBD

In vivo study

TBD

#### 7. Bioactivity

Biological target:

Lavendustin B is an inhibitor of HIV-1 integrase interaction with LEDGF/p75 with an IC50 of 94.07 μM.

In vitro activity

TBD

In vivo activity

TBD

# **Product data sheet**



Note: The information listed here was extracted from literature. MedKoo has not independently retested and confirmed the accuracy of these methods. Customer should use it just for a reference only.