

# Product data sheet



MedKoo Cat#: 572800 Name: NLN 45683 CAS: 2054345-68-3 Chemical Formula: C <sub>14</sub> H <sub>26</sub> N <sub>4</sub> O <sub>7</sub> Exact Mass: 362.1801 Molecular Weight: 362.383		
Product supplied as:		Powder
Purity (by HPLC):		≥ 98%
Shipping conditions		Ambient temperature
Storage conditions:		Powder: -20°C 3 years; 4°C 2 years. In solvent: -80°C 3 months; -20°C 2 weeks.

## 1. Product description:

NLN 45683, also known as N-Boc-Azido-tris(ethylenoxy)-L-alanine is a PEG Linker. This product has no formal name at the moment. For the convenience of communication, a temporary code name was therefore proposed according to MedKoo Chemical Nomenclature (see web page: <https://www.medkoo.com/page/naming>).

## 2. CoA, QC data, SDS, and handling instruction

SDS and handling instruction, CoA with copies of QC data (NMR, HPLC and MS analytical spectra) can be downloaded from the product web page under “QC And Documents” section. Note: copies of analytical spectra may not be available if the product is being supplied by MedKoo partners. Whether the product was made by MedKoo or provided by its partners, the quality is 100% guaranteed.

## 3. Solubility data

Solvent	Max Conc. mg/mL	Max Conc. mM
TBD	TBD	TBD

## 4. Stock solution preparation table:

Concentration / Solvent Volume / Mass	1 mg	5 mg	10 mg
1 mM	2.76 mL	13.80 mL	27.60 mL
5 mM	0.55 mL	2.76 mL	5.52 mL
10 mM	0.28 mL	1.38 mL	2.76 mL
50 mM	0.06 mL	0.28 mL	0.55 mL

## 5. Molarity Calculator, Reconstitution Calculator, Dilution Calculator

Please refer the product web page under section of “Calculator”

## 6. Recommended literature which reported protocols for in vitro and in vivo study

In vitro study

TBD

In vivo study

TBD

## 7. Bioactivity

Biological target:

Azido-PEG3-Ala-Boc is a PEG-based PROTAC linker that can be used in the synthesis of PROTACs.

In vitro activity

TBD

In vivo activity

TBD

*Note: The information listed here was extracted from literature. MedKoo has not independently retested and confirmed the accuracy of these methods. Customer should use it just for a reference only.*